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*Measurements of DNA in microvolumes have become a routine operation in many life science laboratories. As the nitrogenous bases in nucleotides absorb light at 260 nm, this wavelength can be used to determine the DNA concentration. Typical applications are quantification of template prior to sequencing and PCR or measurement of sample purity.*

*Extinction coefficients of nucleic acids are well established:*

*Based on the average extinction at 1 OD of  $0.020 (\mu\text{g/ml})^{-1} \text{cm}^{-1}$  for double-stranded DNA,  $0.027 (\mu\text{g/ml})^{-1} \text{cm}^{-1}$  for single stranded DNA and  $0.025 (\mu\text{g/ml})^{-1} \text{cm}^{-1}$  for RNA, nucleic acid measurements can be done without using a standard curve.*

Berthold Detection Systems GmbH  
Bleichstrasse 56–68  
D-75173 Pforzheim/Germany

Phone: +49(0)7231/9206-0  
Fax: +49(0)7231/9206-50  
E-Mail: [contact@titertek-berthold.com](mailto:contact@titertek-berthold.com)  
Internet: [www.titertek-berthold.com](http://www.titertek-berthold.com)

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*Sample concentration will be automatically calculated by Colibri. Different protocols are available for each nucleic acid type.*

Axit nucleic (Nucleic Acids)	Độ hấp thụ quang riêng (Extinction coefficient) [ $\text{cm}^{-1} \text{M}^{-1}$ ]	Chu trình Colibri (Colibri protocol)
DNA sợi đôi (Double stranded DNA)	50	dsDNA-50
DNA sợi đơn (Single stranded DNA)	33	ssDNA-33
RNA	40	RNA-40

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*Here we report the measurement of ssDNA Oligonucleotide.*

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
*In a first step the concentration of the oligonucleotide stock solution was determined. 2  $\mu\text{l}$  were measured in quadruplicates. Afterwards the stock solution was diluted 1:5 with  $\text{ddH}_2\text{O}$ . A 1:2 dilution series in  $\text{ddH}_2\text{O}$  was prepared. Measurements were taken in triplicates, measurement volume was 2  $\mu\text{l}$  each. Prior to each measurement series the blank was determined by using  $\text{ddH}_2\text{O}$ .*

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
Software settings:


- Protocol: Nucleic Acid, ssDNA-33
- Pathlength: Normal (1mm)

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
*The Colibri microvolume spectrometer shows excellent linearity and sensitivity down to 1 ng/μl ssDNA. Predefined protocols and automatic calculation of sample concentration makes measurement fast and convenient.*

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*Samples and labspace provided by Sarah Becker, KIT Karlsruhe, Germany, Department of Developmental Biology*